



INS LifeGuard

Voice Reach “Talking” Pendant

With Two-Way Voice Communication like a Speaker Phone.



Our **Voice Reach Talking Pendant** is the smallest and least obtrusive in the market. Compact and easy to use, the Talking Pendant summons assistance whenever you need it, at any time of the day or night.

It is designed to allow immediate voice contact during an emergency call. Within the signal range (100-150 meters) of the console unit, the simple press of a button can open up **two-way voice communication** between you and our Emergency Response Centre. The built-in microphone and speaker lets you talk to our nurse right through the pendant like a speaker phone.

Features

- Smallest and least obtrusive talking pendant in the market
- Two-way, hands-free voice communication allows a nurse in our Emergency Response Centre to talk to the user via the built-in microphone and speaker, even if the user isn't near the console unit
- Ultra-low power consumption
- Compatible with the LifeGuard L-401RV unit
- Can be worn around the neck as a pendant

Specifications

- Frequency Available: 1880-1900MHz (Europe) / 1920-1930 MHz (US)
- Max Range: 150 meters in open space
- Power: Lithium 3V battery
- Battery Life: Over 5.77 years
- Operating Temperature: -20° to +50°C
- Humidity: up to 90% non-condensing
- Dimensions: 61 x 42 x 20 mm

Battery life calculated on an average of three activations per month. Minimum order quantity required for Talking Pendant.

Talking pendant is not waterproof. However, it comes with one of our standard transmitters, which is waterproof and can be used in the bath or shower.

Specifications subject to change at discretion of The INS Group.

Warning for those with pacemakers:

Radiofrequency energy from mobile devices and other technologies, including our LifeGuard SmartMobile, SmartWatch and SmartTracker, may cause interference with pacemakers. In particular, note the following precautions:

- **PROPER STORAGE** – avoid placing or wearing them next to the pacemaker. For instance, don't carry your mobile device in your shirt or jacket pocket if it rests near the pacemaker.
- **PROPER DISTANCE** – hold your mobile device to the ear opposite the side of the body where the pacemaker is implanted. Ensure your mobile device remains at least 15 centimetres away from your pacemaker.